# COVID-19 Antigen Schnelltestkit

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## VERWEISE

[1] Lai et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. International Journal of Antimicrobial Agents. 55:3; 2020.

[2] Chen Wang, Peter W Horby, Frederick G Hayden, George F Gao. (2020). A novel coronavirus outbreak of global health concern. The Lancet, 395(10223), 470-473

[3] World Health Organization: Clinical management of severe acute respiratory infection when novel coronavirus (nCoV) is suspected: Interim Guidance. 12 January, 2020.

[4] Clinical and Laboratory Standards Institute. Viral Culture; Approved Guidelines. CLSI document M41-A [ISBN 1562386239] Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne

## SYMBOL INDEX

	For <i>in vitro</i> diagnostic use only		Do not reuse
	Expiry date		See instruction for use
	Warning, please refer to the instructions in the annex		Manufacturer
	Temperature scope within which the product is reserved		Batch number
	Tests per kit		Date of manufacturer
	European union authorized representative		CE mark
	Keep away from sunlight/ Keep away from heat		Keep dry

Beijing Beier Bioengineering Co., Ltd.  
No.99 Chuangxin Road, Lucheng Industrial Development Zone,  
Huangcun Town, Daxing District, Beijing, P. R. China  
Tel: +86 010-61208560, Fax: +86 010-61208569

MedNet EC-REP GmbH  
Borkstrasse 10, 48163 Muenster, Germany

## Swab

GUANGZHOU VICTOR MEDICAL TECHNOLOGY CO.,LTD  
4F, Dalang (Factory A1), TangMei Village, Xintang Town, Zengcheng District, Guangzhou  
511340, China

MedPath GmbH  
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

Revision: January 25, 2021 A02B01

# COVID-19 Antigen Rapid Test Kit

## INTENDED USE

The COVID-19 Antigen Rapid Test Kit is an immunochromatographic intended for the direct and qualitative detection of SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider. The test is to be used as an aid in the diagnosis of coronavirus infection disease (COVID-19), which is caused by SARS-CoV-2.

The COVID-19 Antigen Rapid Test Kit does not differentiate between SARS-CoV and SARS-CoV-2.

Results are for the identification of SARS-CoV-2 nucleocapsid protein antigen. Antigen is generally detectable in upper respiratory specimen during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Negative results do not rule out COVID-19 and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions.

For in vitro diagnostic use only. For professional use only.

## SUMMARY

A novel coronavirus (2019-nCoV) was identified in December 2019, which has resulted in hundreds of thousands of confirmed human infections worldwide. On February 11, 2020 the International Committee for Taxonomy of Viruses (ICTV) renamed the virus SARS-CoV-2. The symptoms of COVID-19 are similar to other viral respiratory diseases and include fever, cough, shortness of breath.

**The SARS-CoV-2 virus has four structural proteins known as S (spike), E (envelope), M (membrane) and N (nucleocapsid) proteins; the N protein holds the RNA genome and the S, E and M proteins together form the viral envelope. Immunological tests currently available on the market make use of the N protein as an antigen, among other things, in order to react with the rapid diagnostic reagents.**

## PRINCIPLE

The COVID-19 Antigen Rapid Test Kit is designed to detect the presence or absence of SARS-CoV or SARS-CoV-2 nucleocapsid proteins by sandwich method. When specimen are processed and added to the sample well, the specimen is absorbed into the device by capillary action. If SARS-CoV or SARS-CoV-2 antigens present in the specimen, it will bind to the SARS-CoV-2 Antibody-labeled conjugated and flows across the coated nitrocellulose membrane in the test strip.

When the SARS-CoV or SARS-CoV-2 antigens level in the specimen is at or above the detection limit of the test, the antigens bound to the SARS-CoV-2 antibody-labeled conjugate are captured by another SARS-CoV-2 antibody immobilized in the Test line (T) of the device, and this produces a colored test band that indicates a positive result. When the SARS-CoV or SARS-CoV-2 antigens level in the specimen does not exist or detection limit of the test, there is not a visible colored band in the Test line (T) of the device. This indicates a negative result.

The nucleocapsid protein is the most conserved antigen in the SARS-CoV-2 structure. To date, no mutants of the nucleocapsid protein have been discovered. Current mutations on the SARS-CoV-2 only occur in the S protein, regardless of whether they are, for example, British or South African variations, these are mutations of the S protein.

## PRECAUTION

1. This kit is for in vitro diagnostic use only.
2. For healthcare professionals and professionals at point of care sites.
3. Do not touch the reaction area of test strip.
4. Please read all the information in this leaflet before performing the test.
5. The test cassette should remain in the sealed pouch until use.
6. All specimen should be considered potentially hazardous and handled in the same manner as an infectious agent.
7. Do not use test kit beyond the expiration date.
8. To avoid erroneous results, specimen must be processed as indicated in the test procedure section.
9. Testing should be applied by professionally trained staff working in certified laboratories or clinics at which the sample(s) is taken by qualified medical personnel. Wear protective clothing such as laboratory coats, disposable gloves, and eye protection.
10. Do not reuse any components of the kit.
11. When collecting a sample, use the swab supplied in the kit. Use of alternative swabs may result in false.
12. The test result should be interpreted by the physician along with clinical findings and other laboratory test results.

13. Disposal of the diagnostic: All specimen and the used-kit has the infectious risk. The used test cassette should be discarded according to federal, state and local regulations.

## MATERIAL

### Material Provided

- 1.20×Sterile Swabs
- 2.20×Extraction tube with Sample Lysis Solution
- 3.20×Test cassette
- 4.1×Instructions for use

### Material Required but Not Provided


- 1.Personal protective equipment, such as protective gloves, medical mask, goggles and lab coat.
- 2.Timer

## STORAGE AND STABILITY


- 1.Store as packaged in the sealed pouch at temperature (4-30°C or 40-86 F).
- 2.Once open the pouch, the test should be used within 30 minutes. Prolonged exposure to hot and humid environment will cause product deterioration.
- 3.The cassette is stable within the expiration date printed on the labeling.
- 4.The LOT and the expiration date were printed on the labeling.

## SPECIMEN COLLECTION AND PREPARATION


### Nasopharyngeal Swab Specimen Collection

1. Tilt patient's head back 70 degrees


2. Insert swab into nostril. (Swab should reach depth equal to distance from nostrils to outer opening of the ear.) Leave swab in place for several seconds to absorb secretions.

3. Slowly remove swab while rotating it. (Swab both nostrils with same swab.)

### Oropharyngeal Swab Specimen Collection

1. For the oropharyngeal swab, take the swab, insert it into the mouth and swab the back of the throat and tonsils. (Avoid the tongue)

### Specimen Transport and Storage

Freshly collected specimen should be processed as soon as possible. The oropharyngeal or nasopharyngeal swabs are stable for up to 24-hours at room temperature or 2° to 8°C.

## TEST PROCEDURE

Please read the instruction for use carefully before performing the test. Allow the Test cassette, Sample Lysis Solution and specimen to equilibrate to temperature (15-30 °C or 59-86°F) before testing.

1. Take an extraction tube from the kit and remove a test cassette from the foil pouch by tearing at the notch. Place them on a level surface.

2. Place the swab sample in the sample lysis solution for 15 seconds and soak. Stir well by turning the swab against the side of the extraction tube 5 times.

3. Remove the swab while squeezing the sides of the tube to extract the liquid from the swab. Discard the swab in biohazard waste.

4. Now screw the extraction tube back on. Then gently turn the extraction tube 5 times.

5. Open the upper cap and transfer the sample with 3 drops (80µl) (ATTENTION) into the sample area of the test strip.

**Note:** If using a frozen sample, the sample must be at room temperature before testing.

**6. Interpret the result after 15-20 min, but no more if 20 minutes are exceeded.**

## RESULT INTERPRETATION

### Positive Result


Colored bands appear at both test line (T) and control line (C). It indicates a positive result for the SARS-CoV or SARS-CoV-2 antigens in the specimen.

### Negative Result

Colored band appear at control line (C) only. It indicates that the concentration of the SARS-CoV or SARS-CoV-2 antigens does not exist or below the detection limit of the test.

### Invalid Result

No visible colored band appear at control line after performing the test. The directions may not have been followed correctly or the test may have deteriorated. It is recommended that the specimen be re-tested.

## QUALITY CONTROL

A procedural control is included in the test. The control line is used as the internal procedural control. The appearance of control line indicate a correct procedure, the absence of control line indicate an inappropriate procedure, specimen volume or degenerative product. Good laboratory practice recommends the use of the control materials. Users should follow the appropriate federal state, and local guidelines concerning the frequency of assaying external quality control materials.

## LIMITATIONS OF PROCEDURE

- The contents of this kit are to be used for the qualitative detection of SARS-CoV or SARS-CoV-2 antigens from nasopharyngeal swab and oropharyngeal swab.
- This reagent is a qualitative assay. It is not designed to determine the quantitative concentration of SARS-CoV or SARS-CoV-2 antigens.
- The accuracy of the test depends on the sample collection process. Improper sample collection, improper sample storage, or repeated freezing and thawing of the sample will affect the test result.
- The results from this test are intended to be an aid in clinical reference only. Each physician must interpret the results in conjunction with the patient's history, physical findings, and other diagnostic procedures.
- Limited by the method of antigen detection reagents, for negative test results, it is recommended to use nucleic acid detection or virus culture identification methods for review and confirmation.
- A false-negative test result may occur if the level of viral antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly; Therefore, a negative test result does not eliminate the possibility of SARS-CoV or SARS-CoV-2 infection.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative test results are not intended to rule in other non-SARS-CoV or SARS-CoV-2 viral or bacterial infections.
- Sensitivity of the test after the first five days of the onset of symptoms has been demonstrated to decrease as compared to a RT-PCR SARS-CoV or SARS-CoV-2 assay. This test detects both viable (live) and non-viable SARS-CoV or SARS-CoV-2. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.

## PERFORMANCE CHARACTERISTICS

### A.Limit of detection (LOD)

The LOD for the COVID-19 Antigen Rapid Test Kit was established using limiting dilutions of a viral sample by Heat inactivation. The material was supplied at a concentration of  $2.6 \times 10^5$  TCID<sub>50</sub>/mL. Using this concentration, the LOD was further refined with a 10-fold dilution series. The last dilution demonstrating 100% positivity was then tested in an additional 20 replicates tested in the same way.

Based on this testing the concentration was confirmed as:  $1.3 \times 10^2$  TCID<sub>50</sub>/mL

### B.Sensitivity and Specificity

474 clinical case samples which include 171 confirmed case samples\* and 303 confirmed excluded case samples\*, were obtained for testing, and then compared the test results between The COVID-19 Antigen Rapid Test Kit's result and the clinical diagnosis. The results of sensitivity and specificity between the two methods are show below.

Reagents	Clinical diagnosis		Total
	Positive	Negative	
Beier Reagents	Positive	1	166
	Negative	302	308
Total	171	303	474

\*Confirmed cases were the patients diagnosed according to the treatment plan and PCR result.

\*Confirmed excluded cases were identified by negative PCR results.

Results analysis: Sensitivity=96.5%(95% CI: 93.7%~99.3%)

Specificity=99.7%(95% CI: 99.0% -100%)

Positive predictive values= 99.4%(95% CI: 98.2% -100%)

Negative predictive values= 98.0%(95% CI: 96.5%-99.6%)

Total consistent: 98.5% (95%CI: 97.4%~99.6%)

### C. Cross-reactivity

Cross-reactivity of the COVID-19 Antigen Rapid Test Kit was evaluated by testing a panel of high prevalence respiratory pathogens that could potentially cross-react with the COVID-19 Antigen Rapid Test Kit. Each organism and virus was tested in triplicate. The final concentration of each organism is documented in the following table.

Cross-Reactivity: COVID-19 Antigen Rapid Test Kit- Wet Testing		
Virus/Bacteria/Parasite*	Concentration	Cross-Reactivity (Yes/No)
Human coronavirus 229E	$1.0 \times 10^5$ U/mL	No
Human coronavirus OC43	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Human coronavirus NL63	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Adenovirus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Human Metapneumovirus	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Parainfluenza virus 1	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Parainfluenza virus 2	$1.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Parainfluenza virus 3	$5.1 \times 10^5$ TCID <sub>50</sub> /mL	No
Parainfluenza virus 4	$1.5 \times 10^4$ TCID <sub>50</sub> /mL	No
Influenza A	$2.7 \times 10^5$ TCID <sub>50</sub> /mL	No
Influenza B	$3.0 \times 10^5$ TCID <sub>50</sub> /mL	No
Enterovirus D68	$3.6 \times 10^5$ TCID <sub>50</sub> /mL	No
Respiratory syncytial virus	$4.2 \times 10^5$ TCID <sub>50</sub> /mL	No
Rhinovirus	$1.0 \times 10^5$ PFU/mL	No
MERS-coronavirus	$1.5 \times 10^5$ TCID <sub>50</sub> /mL	No
Haemophilus influenza	$1.5 \times 10^6$ CFU/mL	No
Streptococcus pneumoniae	$1.1 \times 10^6$ CFU/mL	No

### D.Interferences

The test result of The COVID-19 Antigen Rapid Test Kit do not be interfered with the substance at the following concentration:

Interfering substance	Concentration
Purified Mucin	5%
Human blood	4%
Nasal spray(0.9% NaCl)	150μL
Afrin (Oxymetazoline)	15%
Tobramycin	3mg/dL
Fluticasone	126ng/dL
Budenoside	630ng/dL
Dexamethasone	1.2mg/dL

### E.Precision

Reproducibility studies were performed for the COVID-19 Antigen Rapid Test Kit at three physician office laboratories(POL). One hundred fifty (150) clinical swab specimen including 50 negative, 50 borderline positive and 50 positive, were used in this study. Each specimen was run in triplicate for three days at each POL. The intra-assay agreements were 100%. The inter-site agreement was 100%.

## REFERENCES

- [1] Lai et al. Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) and coronavirus disease-2019 (COVID-19): The epidemic and the challenges. International Journal of Antimicrobial Agents. 55:3; 2020.
- [2] Chen Wang, Peter W Horby, Frederick G Hayden, George F Gao. (2020). A novel coronavirus

outbreak of global health concern. The Lancet, 395(10223), 470-473.

[3] World Health Organization: Clinical management of severe acute respiratory infection when Novel coronavirus (nCoV) infection is suspected: Interim Guidance. 12 January, 2020

[4] Clinical and Laboratory Standards Institute. Viral Culture; Approved Guidelines. CLSI document M41-A [ISBN 1562386239] Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne

## INDEX OF SYMBOL

	For <i>in vitro</i> diagnostic use only		Do not reuse
	Expiry date		See instruction for use
	Warning, please refer to the instructions in the annex		Manufacturer
	Temperature scope within which the product is reserved		Batch number
	Tests per kit		Date of manufacturer
	European union authorized representative		CE mark
	Keep away from sunlight/ Keep away from heat		Keep dry

Beijing Beier Bioengineering Co., Ltd.  
No.99 Chuangxin Road, Lucheng Industrial Development Zone,  
Huangcun Town, Daxing District, Beijing, P. R. China  
Tel: +86 010-61208560, Fax: +86 010-61208569

MedNet EC-REP GmbH  
Borkstrasse 10, 48163 Muenster, Germany

### Swab

GUANGZHOU VICTOR MEDICAL TECHNOLOGY CO.,LTD  
4F, Dalang (Factory A1), TangMei Village, Xintang Town, Zengcheng District, Guangzhou  
511340, China

MedPath GmbH  
Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

Revision: January 25, 2021 A02B01